

# Pig Training — Reference Card

---

Print this page and tape it to the rig.

```
+=====+
|           PIG TRAINING – DAILY USE           |
+=====+
```

STEP 1. HOPPER + MED-PC (Shaping/pellet sessions only)

- 
- [ ] Hopper has ~150+ pellets
  - [ ] MED-PC.exe is running on the rack PC
  - [ ] In MPCLoader, run loader macro -> click S;1
  - [ ] MED-PC screen shows Box 1 with empty counters

STEP 2. LAUNCH

-----

Double-click Start Training on the desktop  
A black terminal window opens AND a small  
"Pig Training" window appears.

STEP 3. PICK THE TASK

-----

```
+-----+
| Task: [ Shaping (MED-PC pellet) v ] |
|           [ Motor Task - Accuracy   ] |
|           [ Motor Task - Reaction Time ] |
|                                     |
|           [ Start Training   ] |
+-----+
```

Read the notes box, then click Start Training.

STEP 4. CONFIGURE

-----

The settings popup opens.

- Pick the stage to start in
- Set responses per stage (defaults are fine)
- Autoshape: Yes (for new subjects)

Click [ Start ].

STEP 5. MONITOR

-----

Watch the terminal window for reward lines:

```
trial started - stage 0
[REWARD] trial 1: requested=1 delivered=1 status=ok
```

Each [REWARD] line = one pellet was delivered.  
Don't touch the keyboard. Don't close anything.

STEP 6. END

- 
- Session completes automatically.
  - Report popup shows stats. Click [ Exit ].
  - Terminal says "Disconnected!" -- press any key to close it.

To stop early: click the tiny gray square in the  
top-left corner of the touchscreen.

STEP 7. DATA

-----

Shaping -> C:\MED-PC\Data\  
Motor Tasks -> HamTraining\Data\\*.csv

Don't move or delete these without the lab tech.

```
+=====+  
|  Stuck?  Don't poke the terminal.  Call the lab tech.  |  
+=====+
```